



POSITIVE POSTURE®

DUALTECH PRO AI

OWNER'S MANUAL

Read this manual thoroughly before you start using your massage chair. In particular, please pay close attention to the sections of this manual devoted to safety.

Follow the instructions in this manual for the safe use and maintenance of your massage chair. This manual should be kept on hand as a reference.

This product is designed for home use.

**POSITIVE POSTURE
DUALTECH PRO AI MASSAGE CHAIR**

CONTENTS

02	Introduction
03	Cautions & Warnings
11	Main Features
13	Parts & Functions
15	Assembly
18	Moving the Chair
19	Sessions
21	Operations & Use
33	Troubleshooting
36	Warranty & After Sales Service
37	Specifications







INTRODUCTION

Welcome to the world of Positive Posture and the transformative experience of your new DualTech Pro AI massage chair. More than just a place to sit, this innovative chair is your personal sanctuary for rejuvenation and well-being. Designed to integrate seamlessly into your life, the DualTech Pro AI offers a pathway to enhanced relaxation, reduced tension, and a renewed sense of balance, allowing you to face each day feeling refreshed and revitalized.

We understand the demands of modern life and have crafted the DualTech Pro AI to be your dedicated partner in self-care. With intuitive features designed to understand your unique needs, the DualTech Pro AI empowers you to prioritize your well-being effortlessly, contributing to a more comfortable, balanced, and ultimately, more enjoyable life. We are confident that your DualTech Pro AI will become an indispensable part of your daily routine, and we sincerely hope you will cherish the exceptional care and lasting benefits it provides.

CAUTIONS & WARNINGS

Indications, Symbols, & Descriptions

 WARNING	This action could result in serious injury or death
 CAUTION	This action could result in personal injury or damage to property
	This action is prohibited
	Further action is required
	Do not disassemble
	Unplug power cord from wall outlet

Important Safety Precautions

- Please follow all basic precautions while using this chair.
- Please read all instructions before using this massage chair. Refer to your Owner's Manual for important information about assembly as well as correct and safe procedures for using the chair.
- The above symbols are used to identify important safety instructions. Please familiarize yourself with these symbols and their meanings.
- Remember to keep this Owner's Manual in an easily accessible place.

CAUTIONS & WARNINGS

General Safety



TO REDUCE THE RISK OF ELECTRIC SHOCK

Always unplug this appliance from the electrical outlet immediately after using and before cleaning.



TO REDUCE RISK OF BURNS, FIRE, ELECTRIC SHOCK, OR INJURY TO PERSONS

- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts.
- Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
- Close supervision is necessary when this appliance is used by, on, or near children or individuals with disabilities.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this appliance if the cord or plug is damaged, if it is not working properly, if it has been dropped or damaged, or if it has been exposed to water. If this happens, contact the service center for examination and repair.
- Do not carry this appliance by the supply cord or use the cord as a handle.
- Keep the cord away from heated surfaces.
- Never operate the appliance with air openings blocked. Keep the air openings free of lint, hair, etc.
- Never drop or insert any object into any opening.
- Do not use outdoors.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- To disconnect, turn all controls to the off position, then remove plug from the outlet.
- Use heated surfaces carefully. May cause serious burns. Do not use over areas with sensitive skin or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
- Keep children away from extended foot support (or other similar parts).
- Connect the appliance to a properly grounded outlet only. See Grounding Instructions.

CAUTIONS & WARNINGS

Grounding Instructions

This product must be grounded. If it should malfunction or break down, grounding provides a path of least resistance for electric shock. This product is equipped with an equipment-grounding conductor and a grounded plug. The plug must be plugged into an appropriate outlet that is properly installed and grounded in accordance with all local codes and ordinances.

DANGER: Improper connection of the equipment-grounding conductor can result in a risk of electric shock. Check with a qualified electrician or serviceman if you are in doubt as to whether the product is properly grounded. Do not modify the plug provided with the product—if it will not fit the outlet, have a proper outlet installed by a qualified electrician.

This product is for use only on a nominal 120V circuit and has a grounded plug that looks like the plug illustrated in Figure A below. Make sure that the product is connected to an outlet having the same configuration as the plug. An adapter should never be used with this product, such as shown in Figure B below.




Figure A




Figure B

CAUTIONS & WARNINGS

Safety During Installation



ELECTRICAL
GROUNDING

- This product uses a normal Alternating Current (AC) circuit and has a ground pin as displayed in the figure to the right. Always make sure that the chair is connected to an outlet that has the same configuration as the plug, and is fully inserted into the outlet. If the chair should malfunction or break down, the grounding will provide a path of least resistance for electrical current to help reduce the risk of electric shock.
- This chair must be plugged into an appropriate outlet that has been properly installed and grounded in accordance with all local codes and ordinances. We strongly recommend using a high-quality surge protector with your chair. A surge protector offering at least 2,000 joules of surge suppression is required for warranty service.
- Do not use with any type of transformer device.
- Do not use the chair with a power supply other than those shown in the image on page 5 (Model, Power Supply AC120V).
- Make sure the main switch is in the OFF position before inserting the plug into the outlet.
- Any modification of the plug provided is prohibited and can cause serious injury. If the plug does not fit into an outlet, contact a qualified electrician for proper installation.
- Failure to connect the grounding conductor of the chair can result in the risk of electrical shock. If you feel the product is not properly grounded, please contact a qualified electrician.
- Please unplug this chair from the electrical outlet after use and before any cleaning or maintenance. To do so, ensure the main power switch is in the OFF position, then unplug the chair from the outlet.



ELECTRICAL
GROUNDING

- This product is intended for use described in this Owner's Manual. Do not use any attachments that are not recommended by the manufacturer.
- Do not use this appliance if the cord or plug is damaged, if it is not working properly, if it has been dropped or damaged, or if it comes into contact with water. In case of any of these scenarios, please contact Positive Posture Service immediately.
- Make sure all air openings are free of lint, hair, or other debris. Do not operate this appliance with the air openings blocked.
- Never use this appliance where aerosol spray products are used or oxygen is being administered.
- Excessive heating of this product can occur and cause fire, electric shock, or serious injury. Do not operate under a blanket, pillow, or heated surface.
- Do not try to move or carry this appliance by pulling or holding onto the power cord.
- Keep power cord away from all heated surfaces.

CAUTIONS & WARNINGS

Safety During Installation



TO AVOID ACCIDENT AND INJURY

- Install the chair on a flat and level surface.
- Allow at least 3" (7 cm) clearance from the chair's back to a wall to allow sufficient space for reclining up and down.
- Do not drop or slam the footrest.
- This appliance is not recommended for children under the age of 14. Please keep children away while in use or storage.

TO AVOID DAMAGE

- Do not use outdoors.
- Stoves, heating products, direct sunlight, or other high temperature areas can cause fading, discoloration, or hardening of the upholstery. Please keep your chair away from this type of exposure.

Safety Before Use



PEOPLE NOT SUITABLE TO USE THE CHAIR

- People with circulatory problems (such as varicose veins or thrombosis), children, and those advised not to use a massage chair by their doctor should not use this chair.
- If you have any of the following health conditions, please consult your doctor for medical advice before using this chair: serious heart problems, serious diabetes, serious osteoporosis, serious skin conditions, malignant tumor(s), sensory impairment, pregnancy (or suspected pregnancy) or those who have just given birth, acute gastrointestinal complaints (gastritis, hepatitis, enteritis), joint dysfunction such as rheumatism, hammer toe and gout, thecitis (or suspected thecitis), high fever, spinal abnormalities due to past injury or ailment, curvature of the spine (scoliosis acute neck sprain (whip lash), or any other medical condition that may contraindicate usage of this product.
- Seek medical advice before using this chair if you have a pacemaker or other electronic medical device, if you have been restricted to bed rest, or if you are planning to use this chair to treat a specific ailment, recover from surgery, or therapy.
- Elderly people and those with weak bones should seek medical advice before using this chair even if they do not have a specific disease or disorder mentioned above.
- Do not use this chair if you have recently consumed alcohol.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.

CAUTIONS & WARNINGS

Safety Before Use, continued...



TO AVOID ACCIDENT, INJURY, AND DAMAGE

- Children should be supervised to ensure that they do not play with the appliance.
- The appliance has a heated surface. Persons sensitive to heat must be careful when using the appliance.
- Do not use this chair for any other reason than described in this manual.
- Do not use this chair with any other therapeutic device or accessory besides those included.
- This chair is intended to be used with all the fabric attached. Failure to do so can result in serious injury.
- Check for any holes or tearing in the active massage area of the chair before each use. If a tear is noticed, stop using immediately and contact Positive Posture Service.
- Check the position of the roller mechanism prior to sitting down. If the rollers are not in the storage position, turn the power on and check that the chair is operating properly. For the proper storage position, see page 15.
- Do not move the chair while someone is in it.
- Do not stand on this chair.
- Do not sit on any part of the chair other than the seat.
- Do not sit, stand, or lean on the footrest unless you are sitting in or getting up from the chair.

Safety During Use



TO AVOID ACCIDENT AND INJURY

- Do not wrap the power or remote cord around your body or place the chair on top of the cord.
- Do not scratch, tear, treat, twist, stretch, or bend the power or remote cord.
- Spilling water or other liquids on this chair may result in electric shock, short circuit, or combustion.
- Keep all children and pets off the chair and away from the surrounding moving area.
- Do not allow children or individuals with disabilities near the chair without close supervision. Keep all persons away from the reclining backrest and extending footrest to avoid injury.
- Do not use the chair while more than one person is sitting in it, or when holding a baby or small child.
- Do not use the chair if any part of your body, including your hands, is wet.
- Do not use the chair unclothed.
- Do not use this chair for any purpose other than its intended uses outlined in this manual.
- Do not use the chair to massage the head, chest, stomach, and/or joints including knees and elbows.

CAUTIONS & WARNINGS

Safety During Use, continued...



TO AVOID ACCIDENT AND INJURY

- Do not place your hands or feet into the roller tracks, or into any moving parts, such as the reclining mechanism.
- Remove all hard objects, such as hair accessories or jewelry, before using.
- Do not use the legrest for anything other than legs or feet.
- Do not operate this massage chair without someone in it.
- Do not put your hands or head between or beneath the legrest.
- Do not sleep in the chair.
- Never unplug the power cord during operation.
- Upon initial use of this chair, start off with a lighter massage. Starting off with a deep massage may result in injury. Elderly and people with weak bones should be careful while selecting massage intensity.
- If you notice any abnormal condition or in case of an emergency, push the Pause button on the remote to stop operation immediately.
- Immediately discontinue use of the chair and seek medical advice if you notice any acute pain or unusual physical discomfort.
- If the massage feels too strong, reduce the stimulation according to this manual. If it still feels too strong, discontinue use immediately.
- Do not use this chair with the back pad removed.



TO AVOID ACCIDENT AND INJURY

- Do not use this chair if there is a power outage or the possibility of one. In this case, stop using the chair immediately, turn the power switch off, and unplug the cord from the wall outlet.
- Do not use this chair if there is a threat of lightning. If so, stop operation immediately, turn the main power switch off, and unplug the chair from the wall outlet.
- If a specific operation does not start or you find any incidence of abnormality, stop using the chair, turn off the main power switch, and unplug the cord from the wall outlet. Call Positive Posture immediately.

CAUTIONS & WARNINGS

Safety During Use, continued...



TO AVOID DAMAGE

- Do not use this chair for any longer than 30 minutes each day. Avoid using a specific focused area massage for more than 5 minutes at a time. Long and continuous massage in the same position may cause discomfort or bodily harm.
- Do not attempt to force your leg, foot, arms, or hands back into the legrest/arm rest if it comes out. Undue force may result in damage to the chair.

OTHER SAFETY PRECAUTIONS

- Do not cause any major impacts to the chair.
- Do not sit on any part of the chair other than the seat.

Safety After Use and When Not in Use



TO AVOID ACCIDENT, INJURY, AND DAMAGE

- Do not plug or unplug the chair with wet hands.
- Do not use any other power cord except the one supplied.
- Grasp the power plug when unplugging the cord from an outlet. Do not use the power cord to pull.
- Turn the main power switch off and unplug the power cord from the wall outlet after each use.
- Always recline the seat back to the upright position after using.
- Wipe down the chair after each use.
- Always store the chair in a low humidity environment.
- Do not place any objects on the chair.
- In case of malfunction, do not attempt to repair. Contact Positive Posture Service.
- Do not modify any part of the chair.
- Be sure to include this Owner's Manual when transferring this chair to another party.
- Do not dispose of this product or its accessories without consulting with your local authorities beforehand.

MAIN FEATURES

WSL-TRACK

An improvement on the SL-Track, Dual Tech Pro AI's WSL-Track ensures consistent, deeply restorative contact from the neck and back to the hips and glutes.

4D² DUAL MECH

Experience the ultimate in precision and personalization with Positive Posture's 4D² Dual Mech technology. Two 4D massage mechanisms work in tandem to deliver simultaneous, full-body therapy. Or adjust them independently so every muscle group gets exactly the care it needs.

INTELLIGENT BODY SCAN

Advanced body scanning technology intelligently maps your unique shape and size at the beginning of each session, locating acupressure points and tailoring each massage session to your body.

AI DE-STRESS TEST™

Positive Posture's proprietary AI Destress Test™ analyzes key wellness indicators—including heart rate, blood oxygen level, and microcirculation—to calculate your personalized Stress Algorithm Score. This advanced algorithm then recommends massage programs best suited to aid relaxation and reduce tension.

20 AUTOMATIC PROGRAMS

This massage chair boasts an impressive selection of 20 automatic programs, offering a diverse range of massage styles and techniques targeting various needs, from deep tissue to gentle relaxation.

FULL BODY AIR MASSAGE

Strategically located air cells provide customizable levels of pressure to improve circulation throughout your body and deliver targeted relief to your shoulders, arms, hands, waist, legs, and feet.

DUAL ZONE HEAT THERAPY

Dual Zone Heat Therapy melts away tension, fatigue, and stiffness as it prepares muscles for massage therapy by applying targeted heat to the arms and lower back.

WINGS OF EMBRACE™ AIR CELLS

Leveraging patented technology, the Wings of Embrace™ shoulder air cells gently extend forward and contour around the shoulders as they inflate, enveloping the upper body in soothing, stabilizing pressure. This precise compression enhances circulation, improves posture, and prepares the upper back for a deeper, more effective massage.

DUAL ACTION REFLEXOLOGY FOOT ROLLERS

Soothe tired feet and stimulate wellness-enhancing pressure points with the touch of two expertly-crafted mechanisms. One travels back and forth along the arch in a firm but gentle arcing motion, while the other rotates specially-shaped massage nodes across the ball of the foot.

DUAL ACTION CALF ROLLERS

Treat your calves to the synchronized touch of two distinct mechanisms. One glides in an up-and-down rhythm to mobilize muscle tissue and improve circulation; the other rotates in a deeply therapeutic circular motion.

THERAKNEAD™ CALF AIR MASSAGE

Soothe your calves with a graceful, rhythmic motion that combines gentle kneading and therapeutic compression. Air cells rotate in a harmonious circular embrace, melting away tension, revitalizing tired muscles, and encouraging healthy circulation—leaving your legs refreshed, renewed, and ready to carry you forward.

MAIN FEATURES

TRUE ZERO GRAVITY® RECLINE

Experience weightless relaxation in True Zero Gravity® recline. This position elevates your legs and evenly distributes your weight, minimizing pressure on your spine and promoting deep relaxation. Float away as you recline into perfect equilibrium, enhancing the benefits of your massage.

HIGH THERAPY POSITION

The high therapy position places your feet well above your heart to improve circulation, reduce swelling and inflammation in the legs, improve respiratory and cardiac function, and hasten injury recovery.

INTUITIVE TOUCHSCREEN CONTROLLER

Effortlessly control your massage experience with our intuitive tablet touchscreen remote. Easily navigate through massage programs, adjust intensity levels, and target specific areas with a simple touch.

QUICK CONTROL DIAL

Turn your DualTech Pro AI on and off, scroll through massage programs, and effortlessly change massage intensity levels with the convenient Quick Control Dial.

BLUETOOTH® SPEAKERS

Relax to your own soundtrack with the chair's built-in Bluetooth® sound system.

HEADREST PILLOW WITH MAGNETIC STORAGE

Designed for personalized comfort, the adjustable headrest pillow perfectly cradles your head, so you can find the perfect position and melt into pure bliss. Magnets hold the headrest in place when not in use.

REMOVABLE BACKPAD

The DualTech Pro AI comes with a convenient removable backpad, providing an additional level of massage intensity customization.

WIRELESS PHONE CHARGER

“Unplug” in more ways than one. DualTech Pro AI's wireless phone charger makes it easy to stay as (dis)connected as you want during your massage.

USB PORT

With DualTech Pro AI's USB charging port, you can recharge your device while you recharge your mind and body.

AMBIENT LIGHTING

Unwind completely and let gentle ambient light enhance your journey to tranquility.


ULTRA SPACE-SAVING DESIGN

Save space in your home with DualTech Pro AI's wall-hugging design, requiring less than 1 inch of space between the chair back and wall.


SLEEK DESIGN

DualTech Pro AI's sleek and stylish design, incorporating wood grain accents, makes it a welcome addition to any living room or den.

PARTS & FUNCTIONS



PARTS & FUNCTIONS



ASSEMBLY

Safety During Installation




TO AVOID DAMAGE


- Do not install the chair in an area of high humidity, such as a bathroom or sauna.
- Do not use the chair in an excessively humid or dusty environment as this may result in malfunction or electric shock.
- Do not use the chair in a room with a temperature of 104-degrees or higher.
- Do not expose the chair to heaters, stoves, or direct sunlight.
- Keep fire sources, such as cigarettes and ashtrays, away from the chair.
- Ensure chair is installed on a flat and level surface to avoid tipping over.
- As the chair is heavy, please exercise extra care not to damage the floor when positioning the chair.

Assembly Preparation


- Ensure there is an allowance of at least 1" (2.5 cm) behind, and 20" (50 cm) in front of the chair to allow the backrest to recline and the legrest to extend without obstruction during operation.
- Maintain a minimum distance of 3' (92 cm) from AV devices, such as a TV or radio, to prevent signal interference.
- It is recommended that a carpet be placed under the chair to prevent scratching the floor and to reduce noise during operation.
- As the chair is heavy, please exercise extra caution.




Base




Leg Rest




Arms (Left & Right)




Remote



Remote Stand



Power Cord



Hardware Bag:
4 Long Screws
2 Short Screws

ASSEMBLY

Assembling the Chair




- To prevent disconnection, be sure to double check all connections by gently pulling on the electrical connectors and air hoses.
- If not connected properly, the air cells will not inflate as designed.

STEP 1

Unbox chair and place in designated area. Carefully remove all items (including footrest from the seat).

STEP 2

Push the remote control cord into the matching connector on the right side near the seat of the chair, then screw until it's tight. Carefully set the remote control on the chair's seat.



STEP 3

Plug the power cord into the back base of the chair (Figure 3A), and turn on the main power switch (Figure 3B). The chair will automatically incline.

Figure 3A






Figure 3B



ASSEMBLY


STEP 4

On the back of the main chair body, find the air lines (one on the left and one on the right) and the heat line (on the left side only). You'll find the matching connectors on the underside of the Back/Seat Pad. Connect the lines by pressing each corresponding line together until they are securely connected. Complete by carefully placing the Back/Seat Pad onto the chair's seat. Press down firmly to secure the velcro.



STEP 5

- Angle the right armrest's four clips into the four holes on the main body of the chair (Figure 5A), also ensuring the metal post on the armrest sits on the hook located next to the shoulder air cell (Figure 5B).
- Hook the bracket (Figure 5C) on the lower part of the armrest over the metal bar located directly below the electrical connectors and airlines on the main body of the chair.
- Push back to slide the arm into place. Make sure there is not a gap between the top of the armrest and the shoulder collar.
- Repeat for left arm.




ASSEMBLY

STEP 6


(Shown with arm detached for visibility)

- Connect the two air lines and the two electrical connectors (*Figure 6A*) from the arm to the main body of the chair. Note: Airlines are color coded, green to green & white to white.
- Install the two longer screws (*Figure 6B*) on the front of the armrest frame.
- Repeat for left arm.



STEP 7


- Line up the leg rest with the front base of the chair, then attach the electrical line and the air line from the leg rest to the chair's base (*Figure 7A*).
- Secure the leg rest by lifting the tops of the metal brackets on the leg rest. Then, carefully lift and place the metal rods on the top of the leg rest onto these brackets (*Figure 7B*). Be gentle to avoid damaging the electrical connection.
- Flip the tops of the metal brackets down over the leg rest and install a short screw into each bracket (*Figure 7C*).



ASSEMBLY


STEP 8

Line up the remote control holder with the slot on the right armrest. Make sure their shapes match. Gently push the holder into the slot. You'll feel and hear a click when it's fully in place. Place the remote control in its holder.



STEP 9

Connect the flap from the top of the legrest to the seat by securing the velcro on the flap on the footrest to the material on the base of the chair beneath the seatpad. Complete by zipping the accordion fabric on the top of the legrest to the seat pad.



STEP 10

Power on the remote control and begin using the chair.


MOVING THE CHAIR



TO AVOID DAMAGE


- Never move the chair while it is operating or when someone is in it.
- As the chair is heavy, please exercise extra care not to abruptly drop it. This may cause damage to the floor as well as the internal components of the chair.

Before moving the chair, turn off the main power and disconnect the power plug from the wall outlet. Bind the cables and put them in the seat to avoid damage when moving the chair.




Option 1

When carrying the chair, the footrest should not be installed. Have one person carefully lift from the front of the chair, while a second person lifts using the handle on the back of the chair. Carry to desired location.



Option 2

To move the chair, first press down firmly on the headrest to shift the weight onto the rear castors. Once the rear castors are engaged, carefully push the chair to your desired location, then gently lower it to the floor.



SESSIONS

Programmed Sessions

Choose from twenty programmed and focused massage sessions. Each session lasts approximately 20 minutes and, unless otherwise noted, can be selected by navigating through the easy-to-use menus on the remote display.

QUICK	Get a taste of some of the unique massage capabilities like 4D ² Dual Mechs, full-body air compression, foot & calf rollers, and TheraKnead™ Calf Air Massage. Quick leaves you feeling centered and renewed and ready to dive deeper into the massage capabilities of the DualTech Pro AI massage chair.
GENTLE	Find abiding serenity in the Gentle massage. This massage session provides a soft, soothing session designed to calm both the body and mind. You'll feel a tranquil release of tension, leaving you deeply relaxed and renewed.
DEEP RELIEF	Unlock your ohhhs and ahhs with the Deep Relief. This massage session is expertly crafted to target and release stubborn knots and tight muscles, delivering intense, focused pressure that goes beyond the surface. You'll experience a profound sense of release as built-up tension melts away, leaving you feeling revitalized and restored.
SPORTS RECOVERY	Optimize your post-workout recovery with Sports Recovery. This massage session is specifically choreographed for the active user, using targeted techniques to help flush out metabolic waste and reduce muscle soreness. It promotes lymphatic drainage and improved circulation, helping your body restore itself faster so you can get back to your routine.
MORNING	Kickstart your day with a Morning massage. This massage session will awaken your body and mind, using well-organized movements that help stimulate circulation and move blood from your extremities back toward the heart. It's the perfect way to get your energy flowing and prepare yourself for the day ahead.
NIGHT	Wind down and slip into sweet slumber with Night. This massage session is designed to prepare your body for sleep, using gentle, repeated movements that help move blood away from your core. It's the perfect way to release the day's tension and drift into a peaceful state of rest.
STRETCH	Realign and mobilize your spine with Stretch. This session is expertly designed to gently pull and extend your body, using a combination of the chair's movements and air cell technology to release tight muscles. You'll experience a lengthening sensation combined with movements to mobilize all levels of your spine that improve flexibility and leaves you feeling limber and open.
NECK & SHOULDER	Release built-up stress with the Neck & Shoulder massage. This massage session is specifically targets the upper back, shoulders, and neck, providing concentrated relief where you hold the most tension. You'll feel a focused release as this area is soothed and revitalized, leaving you feeling lighter and more relaxed.
SPINE ALIGN	Find your balance in Spine Align. This session is carefully choreographed to reduce tension along your spine, allowing your back muscles to relax and your posture to naturally realign. You'll feel a sense of harmony and freedom as stress is released, leaving you feeling centered and restored.
LUMBAR RELIEF	Address low back tension with the Lumbar Relief massage. This massage session focuses attention on the lower back and lumbar region, providing targeted relief to one of the body's most sensitive and pain-prone areas. It works to soothe and release built-up stress, leaving you with a sense of stability and comfort.
FULL BODY 1	Experience a complete refresh with the Full Body 1 massage. This massage session provides a comprehensive, head-to-toe session that methodically soothes every area of your body. You'll feel a total release of tension, leaving you feeling balanced and completely renewed.

SESSIONS

FULL BODY 2	Recharge from head to toe with the dynamic Full Body 2 massage. This program delivers a brisk, all-encompassing session designed to revitalize your entire being. You'll feel an energetic flow as it works every muscle, leaving you feeling awakened and ready for anything.
FULL BODY 3	Indulge in a profound reset with the deep Full Body 3 massage. This massage session thoughtfully works every part of your body with slow, deliberate movements for an all-encompassing, therapeutic experience. You'll feel a deep sense of release, leaving you completely restored and at peace.
LEGS & FEET	Revitalize tired limbs with the Legs & Feet massage. This session provides concentrated relief to your lower body, targeting the legs, calves, and feet. You'll feel an immediate release as the built-up tension from a long day melts away, leaving your lower body feeling energized and restored.
HIPS & GLUTES 1	Relieve lower body tension with the Hips & Glutes 1 massage. This massage session is expertly crafted to target the hips, buttocks, and lower back, providing focused relief to this often-neglected area. You'll feel a deep sense of release as built-up tension melts away, leaving you with a newfound sense of freedom and mobility.
HIPS & GLUTES 2	Unlock deep-seated relief with the targeted Hips & Glutes 2 massage. This massage session delivers slow, deliberate pressure to the lower back and buttocks, working to release stubborn tension in one of the body's most pivotal areas. You'll experience a profound sense of decompression, leaving you feeling balanced and restored from the ground up.
QUIETUDE	Slip into a state of deep calmness with Quietude. This beautifully choreographed session quiets a busy mind and allows your body to gently unwind. You'll feel the weight of your thoughts lift as you sink into quiet relaxation, leaving you with a profound sense of peace and stillness. Yum!
REVITALIZE	Refresh your body and mind with Revitalize. This massage session awakens your sore and tired physique, restoring energy and focus. You'll feel a dynamic release of fatigue, leaving you feeling sharp, alert, and ready for a productive day.
CHINESE MASSAGE	Journey into a timeless tradition with the Chinese massage. This massage session emulates ancient Chinese massage techniques, focusing on balance and flow. You'll feel a harmonious release as the chair works to restore your body's equilibrium, leaving you feeling centered and revitalized.
MY YOGI	Find your flow with the My Yogi. This sweet session is thoughtfully designed to guide your body and mind through a series of slow, deliberate movements, much like a yoga master would. You'll feel a deep sense of presence and release as you move into a state of mindful relaxation and peace.

Manual Sessions

This chair has 7 manual massage methods from which you can choose: Knead, Tap, Pulse, Knead & Pulse, Knead & Tap, Shiatsu, and Roll.

Manual Air Sessions

This chair has 4 manual air sessions from which you can choose: Shoulders & Arms, Waist, Legs & Feet, & Full Body Air.

OPERATION & USE



TO AVOID INJURY

- Always check your surroundings before beginning a massage to ensure that no people, pets, or objects will be hit or squeezed by the reclining back or legrest.



TO AVOID DAMAGE AND INJURY

- This massage chair is meant to be used for indoor household use only. The chair is not designed for commercial purposes.
- If you use the chair in a cold room, do not increase the temperature abruptly. Instead, increase the temperature gradually to a normal level.
- If the chair is stored in a cold place and is moved into a warm environment, wait one hour before use. Otherwise, performance may suffer, or the chair may malfunction due to condensation on the mechanism.
- Do not sit on the chair with the legrest raised. Sudden heavy weight on the legrest could cause damage to the mechanism, or it may result in accident or injury.
- Ensure the backrest returns to the upright position before sitting.
- Do not use the arm or shoulder unit for any area other than your arms and shoulders.
- If your arm comes out of the arm or shoulder unit or your leg/foot come out of the footrest during a massage, do not try to force it back in.
- Do not massage your elbows.
- Do not put anything other than legs and feet in or on the footrest.






TO AVOID DAMAGE

- Do not sit, step, or place heavy objects on the remote control.
- Do not drop the remote control.
- It is recommended that the remote control be placed on the remote stand when not in use.
- Do not exert excessive force on the remote control.
- Do not pull the remote control cord forcefully.
- Do not operate the remote control using sharp edged or hard objects.

OPERATION & USE

Positioning Yourself in the Chair

- Sit back fully in the chair so your entire back is against the backrest and your head is leaning into the pillow.
- Adjust the pillow so it sits just above your shoulders and at the base of your neck. If the pillow is positioned too high or too low, the intelligent body scan will not work correctly. (Figure A)
- Place your arms and elbows between the arm unit air cells with your palms facing down. (Figure B)
- Position your feet in the footrest so your heels and rear calves touch the back of the foot and leg unit. (Figure C)
- Once the leg rest raises, gently press out with your feet to fully extend your legs to a comfortable position. (Figure D)



OPERATION & USE



TO AVOID INJURY

- If you feel the massage is too strong, reduce the stimulation using the Massage Intensity button as described in this manual. If the massage remains too strong, stop operation immediately.
- Do not place hands or feet into the roller tracks, into the space between the seat and back, or into the moving parts.
- Do not sit wearing hard objects, such as hair accessories or jewelry.
- Keep massage sessions to no more than 30 minutes in length. Do not exceed 5 minutes for an area-specific massage for any part of the body.

Retrieving & Storing the Remote

- The remote should be stored on the remote stand when not using the chair.
- It is recommended to leave the remote on the remote stand during a massage, but to retrieve the remote, gently pull it off of the magnetized remote stand.
- To return the remote to its storage position, gently place it on the magnetized stand.

Powering On

- To turn the remote control power on, gently tap the face of the remote's screen OR press and hold the center Power Button located in the center of the dial on the left armrest.
- After powering on, the arm dial will illuminate and the brand logo will briefly display before automatically transitioning to the Programs screen.

Powering Off

- To turn the remote control power OFF, tap the Power button located on the lower left side of the remote screen OR press and hold the Power button located in the center of the dial on the left armrest.
- When the power is OFF, the remote dial is no longer illuminated and the remote screen is black/blank.
- If the chair has been in a reclined position it will automatically restore to the upright position when powered off.

NOTE: The Power button on the remote is not the Main Power switch. See page 13 for Main Power.

OPERATION & USE

Starting a Programmed Massage

- Tap the "PROGRAMS" button on the top of the screen. This will navigate you to the Auto Program menu, where you can select from a variety of pre-set massage programs designed for specific needs and preferences. Note, you can tap the forward/back arrows on the lower left/right of the screen to select from additional program offerings. (See Figure A)
- Tap the desired massage button to activate that massage.
- Once the desired massage is active you will see the active massage summary/options display. (See Figure B)




Figure A




Figure B

Shoulder Scan

To experience an auto massage as intended, it's important for the chair to accurately measure the height of your shoulders prior to the start of a Programmed massage session. To ensure massage accuracy, please follow the steps below: The massage chair will auto-detect your shoulder position by moving from your head down to your shoulder top. During this time, your remote will display a screen indicating this process is in progress. (See Figure A)

Fine Tune Shoulder Position: Once auto-detect is complete, the remote will beep and a screen will appear that allows the roller position to be adjusted by touching and holding the Roller Up or Down arrows on the screen until the mechanism is located in the proper position at the tops of your shoulders. Once positioned correctly, wait 5 seconds for the scan position to register. (See Figure B)

NOTE: Shoulder scan detection does not occur when selecting Manual or Air massages.




Figure A




Figure B

OPERATION & USE

Adjusting Auto Massage Settings

The Auto Program screen appears after a message program has been selected and allows for adjustments to enhance your massage experience.



① Auto Massage Summary Button

Press the Auto Massage Summary button to go back to the summary screen.

② Selected Program Indicator

The name of the currently active message appears here.

③ 4D Intensity Adjustment (Upper Massage Mechanism)

Tap the "Upper Mech Intensity" button to adjust the intensity, or the degree of pressure felt, from the upper massage mechanism. The intensity can be adjusted between five levels, 1 being low and 5 being high. The default setting is 3.

④ 4D Intensity Adjustment (Lower Massage Mechanism)

Tap the "Lower Mech Intensity Adjustment" button to adjust the intensity, or the degree of pressure felt, from the lower massage mechanism. The intensity can be adjusted between five levels, 1 being low and 5 being high. The default setting is 3.

⑤ Air Massage Intensity Adjustment


Tap the "Air Intensity" button to adjust the intensity, or the degree of pressure felt, from the air massage. The intensity can be adjusted between five levels, 1 being low and 5 being high. The default setting is 3.

⑥ Timer

A countdown timer displays the remaining time for your massage.

OPERATION & USE

Starting and Adjusting a Manual Massage



① Manual Button

Tap the “MANUAL” button on the top of the screen. This will navigate you to the Manual Program menu, allowing you to customize your massage experience by manually selecting specific massage techniques, target areas and intensity levels.

② Upper/Lower Mech Selector


Manual massages can be activated and customized for both the Upper and Lower mechanisms. Tap the “UPPER MECH” or “LOWER MECH” button, depending on which mechanism you wish to active/adjust.

NOTE: Both mechanisms have the same massage capabilities and adjustment settings.

OPERATION & USE

③ Mechanism Settings

[Settings appear the same whether upper or lower mech is selected.]



Ⓐ Message Mode

Tap the "Mode" button to activate a message method. Tap repeatedly until the desired Manual method is selected.

Ⓑ Width Adjustment

Width, or the distance between the rollers of the massage mechanism, can be adjusted between narrow, medium and wide by tapping the "Width" button. The default setting is Medium. Feature available for the following Manual methods: Tap, Pulse, Shiatsu, and Roll.

Ⓒ Rolling Range Adjustment

Rolling coverage, or the area of your back that is massaged, can be adjusted between Point or Full by tapping the "Range" button. The default setting is Full. Feature available for all Manual methods.

NOTE: Mechanism will not move when the "Point" range is selected in the "Roll" message mode.

Ⓓ Using the Roller Position Arrows:

In Point, use the up/down "Roller" buttons to move the massage mechanism to the desired location. The up arrow will move the rollers up the track toward the shoulders, while the down arrow will move the rollers down the track toward the glutes.

Ⓔ Speed Adjustment


Speed, or the pace of the massage mechanism movements, can be adjusted in five increments by tapping the "Speed" button. The default setting is 3. Feature available for all Manual methods except Roll & Shiatsu.

Ⓕ 4D Intensity Adjustment

4D Intensity, or the pressure of the massage mechanism, can be controlled in five increments by tapping the "Intensity" button. The default setting is 3. Feature available for all Manual methods.

OPERATION & USE

④ Air Settings



① Air Massage Zone

Tap the "ZONE" button to activate an Air Massage. Tap repeatedly until the desired zone is selected.

② Air Massage Intensity Adjustment

Tap the "Intensity" button to adjust the intensity, or the degree of pressure felt, from the air massage. The intensity can be adjusted between five levels, 1 being low and 5 being high. The default setting is 3.

⑤ Foot Rollers, Calf Rollers, & Heat

① Heat Settings

Tap the "Heat" button to toggle the feature ON or OFF. A highlighted Heat button indicates the feature is active.

NOTE: All Auto programs automatically activate heat.

NOTE: Heat can be enjoyed on its own, or during any massage session.

② Calf Roller & TheraKnead™ Calf Massage Use & Adjustment

Some of the programs automatically activate the calf rollers and TheraKnead™ Calf massage to run intermittently throughout the massage. Tap the "Calf Rollers" button to toggle the feature ON or OFF.

Calf Rollers can also be activated as a manual massage by tapping the "Calf Rollers" button from on the MANUAL mode menu.

NOTE: A highlighted Calf Rollers button indicates the feature is selected and active.


NOTE: Calf Rollers & TheraKnead™ Calf Massage will always activate together, they cannot be turned on and off seperately.

③ Foot Roller Use & Adjustment

Some of the programs automatically activate the foot rollers to run intermittently throughout the massage. Tap the "Foot Rollers" button to toggle the feature ON or OFF.

Foot Rollers can also be activated as a manual massage by tapping the "Foot Rollers" button from on the MANUAL mode menu.

NOTE: A highlighted Foot Rollers button indicates the feature is selected and active.



OPERATION & USE

⑥ Chair Angle Adjustment

During Auto massages, the chair will automatically recline to the suggested position for the massage. If a different reclining angle is desired you can adjust by the following:

Ⓐ Back Position Adjustment

To recline the angle of the chair's back, tap and hold the Back Down button until the desired position is reached. To incline the angle of the chair's back, tap and hold the Back Up button until the desired position is reached. Note adjusting the back position also automatically adjusts the leg position.


Ⓑ Leg Position Adjustment

To recline the angle of the chair's leg, tap and hold the Legs Down button until the desired position is reached. To incline the angle of the chair's leg, tap and hold the Legs Up button until the desired position is reached.

Ⓒ Pre-Set Positions (Zero Gravity & Restore)

The chair is programmed with two pre-set positions that are designed for optimum comfort.

- Zero Gravity: Tap the "Zero Gravity" button to quickly access the Zero Gravity position.
- Restore: Tap "Restore" button to bring the chair to it's full and upright position.



AI De-Stress Test

Using the AI De-Stress Test

- To start, tap the Heart Icon on the bottom left of the remote control, then tap the Start Test button.
- Place your right index finger on the sensor near the front of the right armrest. Insert your finger completely and use light pressure.
- The test takes 30 seconds (a countdown timer will appear on the screen). Remain still until the test is complete.
- After the test, your results will be displayed, along with three recommended massage programs.

NOTE: The screen displays your last five test results. Tap the "Clear Results" button to erase them.



Results Key

- HR: Heart Rate
- SpO2: Blood Oxygen Level
- M-Circ: Micro Circulation
- S-A-S: Stress Algorithm Score

OPERATIONS & USE

Using the Arm Dial & Button


The dial & button located on the left armrest can be used to power the chair, access some of the most popular massage programs, and adjust the 4D massage mechanism intensity without having to navigate through the touch screen remote.

- To Power On: Press and hold the button in the center of the arm dial
- To Power Off: Press and hold the button in the center of the arm dial
- To Activate a Massage: While the chair is powered on, press and quickly release the button in the center of the arm dial. This will activate the Quick Massage. Repeat the quick press and release action to rotate through additional massages (Stretch, Morning & Full Body 3).
- To Adjust 4D Mechanism Intensity: Turn the dial to the right to increase intensity, and to the left to decrease it. This can be used while enjoying auto-programs as well as manual massages.

NOTE: Adjusting the dial will change both the upper and lower mechanism's 4D intensity at that same time.

Pausing a Massage

To pause a massage, tap the "Pause" button, located next to the power button, on the lower left side of the remote screen. Upon doing this the massage mechanism will immediately cease all movement, the chair angle will remain in the exact position it was when the button was pressed. To resume the massage, tap the "Pause" button again (now appearing as a play button). Upon doing this the massage will resume from where it was stopped.





Ending a Massage

To end a massage in the middle of an active session, tap the "Power" button on the lower left side of the screen or press and hold the Power button located in the center of the dial on the left armrest.. This will stop all massaging action, the rollers will move to their storage positions, and the chair will reset to it's fully upright position.



Note: This is the remote power, not the main power switch. Refer to page 23 for the main power switch.



OPERATION & USE

Remote Sound

The remote can be set up to either beep or remain silent with each button press. To change the Sound setting, tap the "SETTINGS" button on the top of the screen. This will navigate you to the Settings menu, where you can select sound ON or Off by tapping the Sound button.

NOTE: The default setting is Sound Off.

LCD Brightness

To adjust the brightness of the LCD display, Tap the "SETTINGS" button on the top of the screen. This will navigate you to the Settings menu, where you can adjust the screen brightness by tapping the LCD Brightness button. Screen brightness can be adjusted between three settings, 1 being dim and 3 being bright. The default setting is 3.

Speakers & Connecting to Bluetooth

- To use the speakers in the massage chair, connect a smart device to the chair's Bluetooth by doing the following:
- Ensure the Main Power is ON.
- From the settings menu of your smart device, select to connect to the Bluetooth connection titled DUALTECHPROAI.
- Your Bluetooth device will indicate once the connection is established.
- Once a smart device is connected, use your smart device to select, play, pause audio and adjust the volume of sound.
- To connect a new device, the currently connected device must first be forgotten which can be done through the settings menu of your smart device.

USB Charging Port

Your chair is equipped with a convenient USB charging port allowing you to easily charge your compatible devices. To use the USB charging port, simply lift the rubber flap that covers the USB port (located on the right armrest) and plug in a compatible cord. Look for the charging indicator on your phone's screen to confirm that charging has started.

NOTE: Charging requires the chair's main power to be on. Make sure the chair is plugged in, and the Main Power Switch is in the ON position.

Accent Lighting (Exterior Armrests)

- Your chair features integrated accent lighting located on the exterior of the armrests. This subtle illumination enhances visibility and adds a touch of modern elegance.
- It is designed for seamless operation and requires no manual control. The default setting is directly linked to the chair's remote power system. When the chair's remote power is activated (i.e. the chair is on and responsive to the remote control), the accent lighting on the armrests will automatically illuminate. When the chair's remote power is turned off (i.e., the chair is off and unresponsive to the remote control), the accent lighting will automatically switch off.

OPERATION & USE

Storage

- After every use, turn off the main power switch located on the back of the chair and unplug the power cord from the wall outlet.
- When not in use, the chair should always be stored in its full upright position (backrest fully inclined, footrest to its lowest position).
- The chair should be stored in a place of low humidity.
- Do not store the chair in direct sunlight or in places of high temperature. This may cause the color of the upholstery to fade.

Cleaning & Care



TO AVOID INJURY

- Before cleaning, always unplug the power cord from the wall outlet.
- Never clean the chair or remote with wet hands.



TO AVOID DAMAGE

- Before cleaning, always unplug the power cord from the wall outlet.
 - Never clean the chair or remote with wet hands.
-
- Use a soft, dry cloth to wipe dust off of the remote control and the polyurethane surfaces located on the inside of the chair.

TROUBLESHOOTING

Troubleshooting Guide

If you encounter problems with the operation of your massage chair, please check the following points to see if they can be easily rectified. If the chair still malfunctions or does not work, please turn off the main power switch and remove the power plug from the wall outlet.

Contact Positive Posture for service. Do not dismantle or attempt to repair the chair yourself.

Doing so will void the warranty.

Problem: My chair doesn't work after powering on.

Possible Cause & Solution:

- Ensure the chair is securely plugged into the wall outlet.
- Ensure the main power switch on the back of the chair is switched to the ON position.
- Ensure the remote control and power cord are fully installed.
- Ensure the remote power is on by tapping the face of the remote screen.
- Ensure a massage session has been selected by touching any programmed or manual massage option on the remote.
- Other issues may include the wires or plug being damaged, the fuse having burned out, or an issue with the inner circuit. Call Positive Posture for service.

Problem: The power cord or plug is abnormally hot.

Possible Cause & Solution:

- Stop using the chair immediately. Turn off the main power switch and remove the power plug from the wall outlet. Call Positive Posture for service.

Problem: My chair is very noisy.

Possible Cause & Solution:

- There are many sounds made by the massage chair that are normal. This includes noise during tapping, the sound of the rollers rubbing against the seat covering while kneading, the sound of the rollers moving up and down the track, and the sound of the air pump and the air bags inflating and deflating. If you hear sounds other than this, or the sounds suddenly become different/louder, please call Positive Posture for service.

TROUBLESHOOTING

Troubleshooting Guide, continued...

Problem: The rollers suddenly stopped in the middle of a session.

Possible Cause & Solution:

- The chair may have lost power. Turn off the main power switch, unplug from the wall, then reinsert the power cord into the wall outlet and turn on the main power switch.
- Ensure the massage session didn't complete its cycle by using the remote to select a new massage session.
- Sometimes when there is too much pressure on the rollers, they can stop moving temporarily. To test this, slightly lift your body or slightly raise the reclining seat back.

Problem: My chair makes a beeping sound at the start of a massage.

Possible Cause & Solution:

- This sound is normal. There is no cause for concern.

Problem: The body scan is not working properly.

Possible Cause & Solution:

- Ensure you are sitting still during the shoulder scan. If you wiggle or move around, you can trigger a false reading.
- Try sitting further back in the chair so that your head rests firmly on the head pillow. You can also adjust the shoulder position with the Roller buttons immediately after the scan is complete.

Problem: The backrest or legrest will not go down entirely.

Possible Cause & Solution:

- Ensure there are not any objects blocking its path.

Problem: The pillow is too high/too low.

Possible Cause & Solution:

- Simply push up or pull down on the pillow to change its position.

Problem: The mechanism is not moving when the manual "Roll" massage mode is active.

Possible Cause & Solution:

- If the "Point" range is selected while the "Roll" massage mode is active, the mechanism will not move.
- Select a different range or massage mode to reactivate the mechanism.

TROUBLESHOOTING

Troubleshooting Guide, continued...

Problem: My chair does not incline when I use the back up button.

Possible Cause & Solution:

- Press and hold the back up button, the chair should return to its upright position within this time.
- Turn off the main power switch and unplug the chair from the wall. Then re-insert the power cord and turn on the main power switch. The chair should return to its upright position.
- If the chair does not, turn the power on with the remote button, then turn the remote power off again. The chair should return to its upright position.
- If the problem persists, please contact Positive Posture service for further assistance.

Q: Can I put the chair on a wood floor?

A: Yes. However, because the chair could damage the floor surface we recommend laying a carpet or something similar under the chair.

Q: I spilled water on the chair. What should I do?

A: There is danger of electric shock. Stop using the chair immediately and do the following:

- Turn off the main power switch on the back left of the chair.
- Unplug the power from the wall outlet.
- Contact Positive Posture for assistance.

Q: How do I make the massage stronger?

A: Try any of the following things:

- Adjust intensity to select maximum pressure.
- Sit back further in the chair so that your weight is against the reclining seat back.
- Lower the reclining seat back.

Q: How do I make the massage softer?

A: Try any of the following things:

- Adjust intensity to select minimum pressure.
- Raise the reclining seat back.

WARRANTY & AFTER SALES SERVICE

Warranty

We're confident that your new massage chair will provide you with many years of relaxation and trouble-free performance. Your Positive Posture massage chair includes a 3-year limited warranty: the first year covers on-site service (in most US/Canada locations) and moving parts, while non-moving parts and the frame are covered for the full three years. Some restrictions apply. See www.PositivePosture.com for full warranty details. We recommend that you print the limited warranty information at the time of purchase and keep it with your Owner's Manual.

Service & Repairs

Before calling for repairs or service, please consult the Troubleshooting Guide. If the problem persists, contact Positive Posture for assistance. For fastest service, please have the following information ready:

- Your name, email, address & telephone number
- Model Number (see the specification field on the front cover of this Owner's Manual or review the sticker affixed to the back base of your chair)
- Serial Number (see the sticker affixed to the back base of your chair)
- Copy of original receipt and purchase date
- Malfunction or issue with as much detail as possible

SPECIFICATIONS

Model Name	DualTech Pro AI
Model Number	PP7605
Power Supply	AC120V 60Hz
Power Consumption	300W
Rated Time	30 minutes
Back Stroke Range	36" (91.5 cm)
Preset Sessions	20 sessions
Manual Sessions	7 session on both upper and lower mech, plus Foot Rollers, Calf Roller & Full Body Air
Air Pressure	Approximately 32 kPa
Program Timer	Approximately 20 minutes
Footrest Adjustment Range	Approximately 7.5" (19 cm)
Reclining Angle	Approximately 135° to 160° against flat floor
Chair Dimensions	Upright: 63.7" (162cm) L x 32.3" (82cm) W x 42.25 (120cm) Reclined: 72" (183cm) L x 32.3" (82cm) W x 35.5" (90cm) H
Seat Width	20.75" (53cm)
Seat Depth	20.75" (53cm)
Chair Weight	293.5 lbs (135.2 kg)
Supplied Items	Power cord, Remote Stand, Remote Control, Leg Rest
Shipping Dimensions	Box 1: 57.9" (147cm) L x 33.1" (84cm) W X 38.2" (97cm) H Box 2: 51.8" (131.5cm) L x 16.6" (42cm) W x 24" (61cm) H
Shipping Weight	Box 1: 315lbs (143 kg) Box 2: 41 lbs (18.5 kg)
Recommended User Weight	265 lbs (120 kg) (maximum)
Recommended User Height	5' 0" (153 cm) to 6' 2" (188 cm)

DUALTECH PRO AI

by  POSITIVE POSTURE®

Your massage chair was meticulously designed to conform to your specific comfort and relaxation needs. If you have questions about the operation of your massage chair, contact information for a Positive Posture chair expert can be found at www.PositivePosture.com.

Positive Posture. A division of Furniture For Life.

© 2025 Positive Posture
Boulder, Colorado | USA
www.PositivePosture.com